



Fully polished RW Crystal Grey granite balance benchtop

Crystal Grey granite features a mid-tone grey background and flecks that add interest and allure to any laboratory space. With the average density between 2.65 and 2.75 g/cm³ The weight of 80mm thick granite is substantial making it a preferred choice for Westlab Spaces weighing tables.

The Westlab Spaces Laboratory Weighing Table Granite is also hard and durable enough to resist abrasion and bear significant weight. It is also considerably inert in nature and can be polished to achieve a brilliant finish

Specifications

Material	Granite
Thickness	80mm
Finishes	Polished
Colour	Crystal Grey
Size (Length and Width)	Cut to size
Compressive Strength	200 MPa [^]
Melting Point (Dry Granite)	1215–1260 °C (2219–2300 °F)
Average density	2.65 and 2.75 g/cm ³
Other Names	Solid Granite Benchtop, Weighing Table, Balance Bench

COLOUR



Crystal Grey

MORE INFORMATION

Order online [here](#)

